



CHIRAL APPLICATION NOTE

Sample name: Dichloroprop

Chromatogram:



Structure:



Run Conditions:

Column: TeicoShell (10cmx4.6mm, 2.7 μ m)
Mobile phase: 100/0.3 wt %: Methanol/Ammonium Trifluoroacetate
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 μ L
Detection: UV 220 nm
Retention (min): 1.22/1.54 (Peak1/Peak 2)

Resolution: 4.91