



CHIRAL APPLICATION NOTE

Sample name: 2-chloronicotine

Chromatogram:



Structure:



Run Conditions:

Column: VancoShell (10cmx4.6mm, 2.7 μ m) x2
Mobile phase: 90/10: Methanol/16 mM Ammonium Formate (pH 3.6)
Flow Rate: 0.7 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 μ L
Detection: UV 263 nm
Retention (min): 6.30/6.57 (Peak 1/Peak 2)

Resolution: 1.57