



CHIRAL APPLICATION NOTE

Sample name: Mandipropamid

Chromatogram:



Structure:



Run Conditions:

Column: TeicoShell (10cmx4.6mm, 2.7 μ m)
Mobile phase: 30/70: Methanol/50 mM Ammonium Formate (pH 3.6)
Flow Rate: 0.5 mL/min
Temperature: 45 °C
Injection Volume: 1.0 μ L
Detection: UV 230 nm
Retention (min): 7.65/8.22 (Peak 1/Peak 2)

Resolution: 1.78